

Setting up and running the Zebrafish Segmentation

Table of Contents

Installation

- install Anaconda (Python) and Libraries
- Setting Up Docker and segmentation model

Run the Segmentation Workflow

- Running the `sort_images_nd2.py` Python Script
- Running the `segmentation.sh` Shell Script to run the Segmentation Model in a Docker Container

Scripts

sort_images_nd2.py: this script is designed to convert Nikon NIS-Elements ND2 files into individual TIFF files. It processes each series, timepoint, channel, and Z-stack within the ND2 file and saves them as TIFF images. The output is organized into directories for each series, with separate subdirectories for bright-field ($C = 0$) and fluorescence images ($C > 0$).

segmentation.sh: Starts a Docker container to run the segmentation model. A path to the folder containing the 'Images' and 'Fluor' folders, created with the previous script, needs to be provided.

[!note] This guide is written for the setup on Mac or Linux. For the use on Windows the `segmentaion.sh` script needs to be exchanged with the `segmentation.ps1`. Other steps can also be slightly different from the description here.

Installation

Install Anaconda and Libraries

Step 1: Install Anaconda

1. **Download Anaconda:**
 - Visit the [Anaconda website](#) and download the Anaconda Installer for macOS.
2. **Install Anaconda:**
 - Open the downloaded `.pkg` file and follow the on-screen instructions to install Anaconda.
3. **Verify Installation:**
 - Open Terminal and type `conda --version` to verify Anaconda is installed.

Step 2: Create a Conda Environment (Optional)

It's a good practice to create a new environment for your projects. To create one, run:

```
conda create --name myenv python=3.8
```

Replace `myenv` with your desired environment name.

Activate your new environment by running:

```
conda activate myenv
```

Step 3: Install `nd2reader` and `skimage` Libraries

With an environment activated (or in the base environment), install the required libraries using `conda`:

```
conda install -c conda-forge nd2reader scikit-image
```

This command installs `nd2reader` and `scikit-image` (which includes `skimage`) from the Conda-Forge repository.

Setting Up Docker

Step 1: Download Docker Desktop

1. Go to the [Docker Hub](#) and download the Docker Desktop application for Mac.
2. Open the downloaded `.dmg` file and drag the Docker icon to your Applications folder to install.

Step 2: Launch Docker

Open Docker from your Applications folder. The first launch might take a while as it sets up. Once Docker is running, open a Terminal and type:

```
docker --version
```

This command will display the Docker version, confirming it's installed.

Step 4: Get the segmentation model as a docker image

The following command downloads and installs our segmentation model:

```
docker pull branhongweili/dqbm_cell_seg:v3.1
```

Run the Segmentation Workflow

Running the `sort_images_nd2.py` Python Script

Open Terminal and activate your python environment by running:

```
conda activate myenv
```

Navigate to the directory containing the `sort_images_nd2.py` script:

```
cd path/to/the/folder/containing/the/script
```

Run the following with the path to your `.nd2`-file:

```
python3 sort_images_nd2.py path/to/your/nd2file
```

The script creates a position folder for each imaged fish, with a folder named 'Images' containing all the bright field images and a folder named 'Fluor' containing all the fluorescence images.

Running the `segmentation.sh` Shell Script to run the Segmentation Model in a Docker Container

Make the Script Executable (only once): In Terminal, navigate to the directory containing the script and run: `bash chmod +x segmentation.sh`

Run the Script: Execute the script by running the following and providing the path to one of the imaged positions (created in the previous step): `bash ./segmentation.sh path/to/position/folder`

This command starts a Docker container, runs the segmentation and saves the resulting segmentation in the position folder.